




Genome polyprotein Antibody

Product Code	CSB-PA189574ZA01Ba
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99DQ9
Immunogen	Recombinant Japanese encephalitis virus Genome polyprotein,partial (2780-3428aa)
Species Reactivity	Japanese encephalitis virus
Tested Applications	ELISA, WB
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein G
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Japanese encephalitis virus
Target Names	Genome polyprotein

Image



Western Blot

Positive WB detected in Recombinant protein
All lanes: Genome polyprotein Antibody antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 87.3 kDa

Observed band size: 105 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.