



Factor Xa inhibitor BuXI Antibody

Product Code	CSB-PA307110ZA01BFM
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P83594
Immunogen	Recombinant Bauhinia unguolata Factor Xa inhibitor BuXI protein
Raised In	Rabbit
Species Reactivity	Bauhinia unguolata
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Factor Xa inhibitor BuXI
Immunogen Species	Bauhinia unguolata (Orchid tree)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	BuXI
Product Modify	Unmodified