



Gellan lyase Antibody

Product Code	CSB-PA308328ZA01GFM
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P85513
Immunogen	Recombinant Geobacillus stearothermophilus Gellan lyase protein
Raised In	Rabbit
Species Reactivity	Geobacillus stearothermophilus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Gellan lyase (EC 4.2.2.25) (Fragments)
Immunogen Species	Geobacillus stearothermophilus (Bacillus stearothermophilus)
Lead Time	Made-to-order (12-14 weeks)
Product Modify	Unmodified