



gltX1 Antibody

Product Code	CSB-PA308728ZA01HUV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P96551
Immunogen	Recombinant Helicobacter pylori gltX1 protein
Raised In	Rabbit
Species Reactivity	Helicobacter pylori
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glutamate--tRNA ligase 1 (EC 6.1.1.17) (Glutamyl-tRNA synthetase 1) (GluRS 1), gltX1
Immunogen Species	Helicobacter pylori (strain ATCC 700392 / 26695) (Campylobacter pylori)
Lead Time	Made-to-order (12-14 weeks)
Target Names	gltX1
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.