



Superoxide dismutase [Mn], mitochondrial Antibody

Product Code	CSB-PA309966ZA01AZV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P86254
Immunogen	Recombinant <i>Alternaria alternata</i> Superoxide dismutase [Mn], mitochondrial protein
Raised In	Rabbit
Species Reactivity	<i>Alternaria alternata</i>
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Superoxide dismutase [Mn], mitochondrial (EC 1.15.1.1) (Manganese-dependent superoxide dismutase) (MnSOD) (allergen Alt a MnSOD) (Fragment)
Immunogen Species	<i>Alternaria alternata</i> (<i>Alternaria rot fungus</i>) (<i>Torula alternata</i>)
Lead Time	Made-to-order (12-14 weeks)
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.