



VTG1 Antibody

Product Code	CSB-PA311247ZA01CH
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P87498
Immunogen	Recombinant Gallus gallus VTG1 protein
Raised In	Rabbit
Species Reactivity	Gallus gallus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Vitellogenin-1 (Minor vitellogenin) (Vitellogenin I) [Cleaved into: Lipovitellin-1 (Lipovitellin I) (LVI); Phosvitin (PV); Lipovitellin-2 (Lipovitellin II) (LVII); YGP42], VTG1, VTGI
Immunogen Species	Gallus gallus (Chicken)
Lead Time	Made-to-order (12-14 weeks)
Target Names	VTG1
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.