




degP Antibody

Product Code	CSB-PA314631ZA01ENV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0C0V0
Immunogen	Recombinant Escherichia coli Periplasmic serine endoprotease DegP (degP) (27-474AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA, WB
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein A/G
Isotype	IgG
Clonality	Polyclonal
Alias	Periplasmic serine endoprotease DegP (EC 3.4.21.107) (Heat shock protein DegP) (Protease Do), degP, htrA ptd
Immunogen Species	Escherichia coli (strain K12)
Target Names	degP

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: degP antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 62.8 kDa

Observed band size: 65 kDa

Product Modify Unmodified

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.