



Con-ikot-ikot Antibody

Product Code	CSB-PA316701ZA01DWI
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0CB20
Immunogen	Recombinant Conus striatus Con-ikot-ikot protein
Raised In	Rabbit
Species Reactivity	Conus striatus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Con-ikot-ikot
Immunogen Species	Conus striatus (Striated cone)
Lead Time	Made-to-order (12-14 weeks)
Product Modify	Unmodified