



Major egg antigen Antibody

Product Code	CSB-PA319453ZA01SXQ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P12812
Immunogen	Recombinant Schistosoma mansoni Major egg antigen protein
Raised In	Rabbit
Species Reactivity	Schistosoma mansoni
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Major egg antigen (p40)
Immunogen Species	Schistosoma mansoni (Blood fluke)
Lead Time	Made-to-order (12-14 weeks)
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.