



Pectate lyase 1 Antibody

Product Code	CSB-PA322455ZA01DYO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P18632
Immunogen	Recombinant Cryptomeria japonica Pectate lyase 1 protein
Raised In	Rabbit
Species Reactivity	Cryptomeria japonica
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pectate lyase 1 (EC 4.2.2.2) (Allergen Cry j I) (Major allergen Cry j 1) (Sugi basic protein) (SBP) (allergen Cry j 1)
Immunogen Species	Cryptomeria japonica (Japanese cedar) (Cupressus japonica)
Lead Time	Made-to-order (12-14 weeks)
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.