



NP Antibody

Product Code	CSB-PA323576ZA01ZAB
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P18272
Immunogen	Recombinant Zaire ebolavirus NP protein
Raised In	Rabbit
Species Reactivity	Zaire ebolavirus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nucleoprotein (Nucleocapsid protein) (Protein N), NP
Immunogen Species	Zaire ebolavirus (strain Mayinga-76) (ZEBOV) (Zaire Ebola virus)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	NP
Product Modify	Unmodified