



nfsA Antibody

Product Code	CSB-PA323751ZA01ENV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17117
Immunogen	Recombinant Escherichia coli nfsA protein
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Oxygen-insensitive NADPH nitroreductase (EC 1.-.-.-) (Modulator of drug activity A), nfsA, mda18 mdaA ybjB
Immunogen Species	Escherichia coli (strain K12)
Lead Time	Made-to-order (12-14 weeks)
Target Names	nfsA
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.