



# nfsA Antibody

|                            |  |
|----------------------------|--|
| <b>Product Code</b>        | CSB-PA323751ZA01ENV  |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.  |
| <b>Uniprot No.</b>         | P17117   |
| <b>Immunogen</b>           | Recombinant Escherichia coli nfsA protein  |
| <b>Raised In</b>           | Rabbit   |
| <b>Species Reactivity</b>  | Escherichia coli   |
| <b>Tested Applications</b> | ELISA, WB (ensure identification of antigen)   |
| <b>Form</b>                | Liquid   |
| <b>Conjugate</b>           | Non-conjugated   |
| <b>Storage Buffer</b>      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, pH 7.4                           |
| <b>Purification Method</b> | Antigen Affinity Purified  |
| <b>Isotype</b>             | IgG  |
| <b>Clonality</b>           | Polyclonal   |
| <b>Alias</b>               | Oxygen-insensitive NADPH nitroreductase (EC 1.-.-.-) (Modulator of drug activity A), nfsA, mda18 mdaA ybjB |
| <b>Immunogen Species</b>   | Escherichia coli (strain K12)  |
| <b>Lead Time</b>           | Made-to-order (12-14 weeks)  |
| <b>Gene Names</b>          | nfsA   |
| <b>Product Modify</b>      | Unmodified   |