



2S seed storage albumin Antibody

Product Code	CSB-PA324829ZA01GGV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P19594
Immunogen	Recombinant Glycine max 2S seed storage albumin protein
Raised In	Rabbit
Species Reactivity	Glycine max
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	2S albumin (GM2S-1) (Napin-type 2S albumin 3) [Cleaved into: 2S albumin small chain (Aspartic acid-rich peptide) (Lunasin); 2S albumin large chain (8 kDa methionine-rich protein) (8 kDa MRP)]
Immunogen Species	Glycine max (Soybean) (Glycine hispida)
Lead Time	Made-to-order (12-14 weeks)
Target Names	GM2S-1
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.