



Ice-structuring protein lambda OP-3 Antibody

Product Code	CSB-PA325687ZA01MAQ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P19606
Immunogen	Recombinant Zoarces americanus Ice-structuring protein lambda OP-3 protein
Raised In	Rabbit
Species Reactivity	Zoarces americanus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ice-structuring protein lambda OP-3 (ISP lambda OP-3) (Antifreeze protein lambda OP-3)
Immunogen Species	Zoarces americanus (Ocean pout) (Macrozoarces americanus)
Lead Time	Made-to-order (12-14 weeks)
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.