



Ubc6 Antibody

Product Code	CSB-PA328517ZA01DLU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P25153
Immunogen	Recombinant Drosophila melanogaster Ubc6 protein
Raised In	Rabbit
Species Reactivity	Drosophila melanogaster
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ubiquitin-conjugating enzyme E2-17 kDa (EC 2.3.2.23) (E2 ubiquitin-conjugating enzyme 6) (Ubiquitin carrier protein) (Ubiquitin-protein ligase), Ubc6, Dhr6 UbcD6
Immunogen Species	Drosophila melanogaster (Fruit fly)
Lead Time	Made-to-order (12-14 weeks)
Target Names	Ubc6
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.