



munIM Antibody

Product Code	CSB-PA331901ZA01MLZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P43641
Immunogen	Recombinant Mycoplasma sp. munIM protein
Raised In	Rabbit
Species Reactivity	Mycoplasma sp.
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Modification methylase MunI (EC 2.1.1.72) (Adenine-specific methyltransferase MunI) (M.MunI), munIM
Immunogen Species	Mycoplasma sp.
Lead Time	Made-to-order (12-14 weeks)
Gene Names	munIM
Product Modify	Unmodified