



soxA Antibody

Product Code	CSB-PA332906ZA01BCG
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P23342
Immunogen	Recombinant Bacillus sp. soxA protein
Raised In	Rabbit
Species Reactivity	Bacillus sp.
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Monomeric sarcosine oxidase (MSOX) (EC 1.5.3.1), soxA
Immunogen Species	Bacillus sp. (strain NS-129)
Lead Time	Made-to-order (12-14 weeks)
Target Names	soxA
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.