



ompD Antibody

Product Code	CSB-PA334291ZA01SXB
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P37592
Immunogen	Recombinant Salmonella typhimurium ompD protein
Raised In	Rabbit
Species Reactivity	Salmonella typhimurium
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Outer membrane porin protein OmpD, ompD, nmpC
Immunogen Species	Salmonella typhimurium (strain LT2 / SGSC1412 / ATCC 700720)
Lead Time	Made-to-order (12-14 weeks)
Target Names	ompD
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.