



Major pollen allergen Car b 1 isoform 2 Antibody

Product Code	CSB-PA336516ZA01CEC
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P38950
Immunogen	Recombinant Carpinus betulus Major pollen allergen Car b 1 isoform 2 protein
Raised In	Rabbit
Species Reactivity	Carpinus betulus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Major pollen allergen Car b 1 isoform 2 (Allergen Car b I) (allergen Car b 1)
Immunogen Species	Carpinus betulus (European hornbeam) (Carpinus caucasica)
Lead Time	Made-to-order (12-14 weeks)
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.