




prgl Antibody

Product Code	CSB-PA336905ZA01SXB
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P41784
Immunogen	Recombinant Salmonella typhimurium Protein prgl (prgl) (1-80AA)
Raised In	Rabbit
Species Reactivity	Salmonella typhimurium
Tested Applications	ELISA, WB
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein A/G
Isotype	IgG
Clonality	Polyclonal
Alias	Protein Prgl, prgl
Immunogen Species	Salmonella typhimurium (strain LT2 / SGSC1412 / ATCC 700720)
Target Names	prgl

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: prgl antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 24.9 kDa

Observed band size: 30 kDa

Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.