



GFP Antibody

Product Code	CSB-PA337004ZA01ADO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P42212
Immunogen	Recombinant Aequorea victoria GFP protein
Raised In	Rabbit
Species Reactivity	Aequorea victoria
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Green fluorescent protein, GFP
Immunogen Species	Aequorea victoria (Jellyfish)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	GFP
Product Modify	Unmodified