



MEL1 Antibody

Product Code	CSB-PA356371ZA01SAC
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04824
Immunogen	Recombinant Saccharomyces cerevisiae MEL1 protein
Raised In	Rabbit
Species Reactivity	Saccharomyces cerevisiae
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Alpha-galactosidase 1 (EC 3.2.1.22) (Alpha-D-galactoside galactohydrolase 1) (Melibiase 1), MEL1
Immunogen Species	Saccharomyces cerevisiae (strain ATCC 204508 / S288c) (Baker's yeast)
Lead Time	Made-to-order (12-14 weeks)
Target Names	MEL1
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.