



yciM Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA359577ZA01ENV |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | A0A0E1LWJ7 |
| Immunogen | Recombinant Escherichia coli yciM protein |
| Raised In | Rabbit |
| Species Reactivity | Escherichia coli |
| Tested Applications | ELISA, WB (ensure identification of antigen) |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | Antigen Affinity Purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Lipopolysaccharide assembly protein B (yciM) |
| Immunogen Species | Escherichia coli (strain K12) |
| Lead Time | Made-to-order (12-14 weeks) |
| Target Names | yciM |
| Product Modify | Unmodified |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |