



esxH Antibody

Product Code	CSB-PA363659ZA01MVZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P9WVK2
Immunogen	Recombinant Mycobacterium tuberculosis esxH protein
Raised In	Rabbit
Species Reactivity	Mycobacterium tuberculosis
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	ESAT-6-like protein EsxH (10 kDa antigen CFP7) (CFP-7) (Low molecular weight protein antigen 7) (Protein TB10.4), esxH, cfp7
Immunogen Species	Mycobacterium tuberculosis (strain CDC 1551 / Oshkosh)
Lead Time	Made-to-order (12-14 weeks)
Target Names	esxH
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.