



aglB Antibody

Product Code	CSB-PA381525ZA01AVE
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A2QEJ9
Immunogen	Recombinant Aspergillus niger (strain CBS 513.88 / FGSC A1513) aglB protein
Raised In	Rabbit
Species Reactivity	Aspergillus niger (strain CBS 513.88 / FGSC A1513)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Probable alpha-galactosidase B (EC 3.2.1.22) (Melibiase B), aglB
Immunogen Species	Aspergillus niger (strain CBS 513.88 / FGSC A1513)
Lead Time	Made-to-order (12-14 weeks)
Target Names	aglB
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.