




omp38 Antibody

Product Code	CSB-PA387720ZA01AUD
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A3M8K2
Immunogen	Recombinant Acinetobacter baumannii Outer membrane protein Omp38 (omp38) (20-356aa)
Species Reactivity	Acinetobacter baumannii
Tested Applications	ELISA, WB
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein G
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Acinetobacter baumannii
Target Names	omp38

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: omp38 antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 52.5 kDa

Observed band size: 66 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.