



GL50803_103887 Antibody

Product Code	CSB-PA419783ZA01GFY
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A8BQJ3
Immunogen	Recombinant Giardia intestinalis GL50803_103887 protein
Raised In	Rabbit
Species Reactivity	Giardia intestinalis
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Endoribonuclease Dicer-like (EC 3.1.26.-)
Immunogen Species	Giardia intestinalis (strain ATCC 50803 / WB clone C6) (Giardia lamblia)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	GL50803_103887
Product Modify	Unmodified