



# Non-structural glycoprotein 4 Antibody

<b>Product Code</b>	CSB-PA465285ZA01ROJ
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	A2T3Q0
<b>Immunogen</b>	Recombinant Rotavirus A Non-structural glycoprotein 4 protein
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Rotavirus A
<b>Tested Applications</b>	ELISA, WB (ensure identification of antigen)
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Non-structural glycoprotein 4 (NSP4) (NCVP5) (NS28)
<b>Immunogen Species</b>	Rotavirus A (isolate Monkey/South Africa/SA11-H96/1958 G3-P5B[2]-I2-R2-C5-M5-A5-N5-T5-E2-H5) (RV-A) (Simian Agent 11 (isolate SI/South Africa/H96/58))
<b>Lead Time</b>	Made-to-order (12-14 weeks)
<b>Target Names</b>	NSP4
<b>Product Modify</b>	Unmodified
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.