



nosM Antibody

Product Code	CSB-PA509748ZA01FNT
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	C6FX52
Immunogen	Recombinant Streptomyces actuosus nosM protein
Raised In	Rabbit
Species Reactivity	Streptomyces actuosus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nosiheptide precursor [Cleaved into: Nosiheptide (NOS) (Antibiotic 9671-RP) (Multhiomycin)], nosM
Immunogen Species	Streptomyces actuosus
Lead Time	Made-to-order (12-14 weeks)
Target Names	nosM
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.