



bgIIIM Antibody

Product Code	CSB-PA670450ZA01BRI
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q45489
Immunogen	Recombinant Bacillus subtilis bgIIIM protein
Raised In	Rabbit
Species Reactivity	Bacillus subtilis
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Modification methylase BgIII (M.BgIII) (EC 2.1.1.113) (N(4)- cytosine-specific methyltransferase BgIII), bgIIIM
Immunogen Species	Bacillus subtilis (strain 168)
Lead Time	Made-to-order (12-14 weeks)
Target Names	bgIIIM
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.