



choA Antibody

Product Code	CSB-PA704121ZA01SPX
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P12676
Immunogen	Recombinant Streptomyces sp. choA protein
Raised In	Rabbit
Species Reactivity	Streptomyces sp.
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cholesterol oxidase (CHOD) (EC 1.1.3.6) (Cholesterol isomerase) (EC 5.3.3.1), choA
Immunogen Species	Streptomyces sp. (strain SA-COO)
Lead Time	Made-to-order (12-14 weeks)
Target Names	choA
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.