



VP30 Antibody

Product Code	CSB-PA755264ZA01ZAT
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q77DJ5
Immunogen	Recombinant Zaire ebolavirus VP30 protein
Raised In	Rabbit
Species Reactivity	Zaire ebolavirus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Minor nucleoprotein VP30 (Transcription activator VP30), VP30
Immunogen Species	Zaire ebolavirus (strain Kikwit-95) (ZEBOV) (Zaire Ebola virus)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	VP30
Product Modify	Unmodified