



Conglutin-7 Antibody

Product Code	CSB-PA757666ZA01ANE
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6PSU2
Immunogen	Recombinant Arachis hypogaea Conglutin-7 protein
Raised In	Rabbit
Species Reactivity	Arachis hypogaea
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Conglutin-7 (2S protein 1) (Seed storage protein SSP1) (Seed storage protein SSP2) (allergen Ara h 2)
Immunogen Species	Arachis hypogaea (Peanut)
Lead Time	Made-to-order (12-14 weeks)
Product Modify	Unmodified