



FXN Antibody

Product Code	CSB-PA810224ZA01MOV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8HXX9
Immunogen	Recombinant Macaca fascicularis FXN protein
Raised In	Rabbit
Species Reactivity	Macaca fascicularis
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Frataxin, mitochondrial (Fxn) (EC 1.16.3.1) [Cleaved into: Frataxin intermediate form; Frataxin mature form], FXN, FRDA1
Immunogen Species	Macaca fascicularis (Crab-eating macaque) (Cynomolgus monkey)
Lead Time	Made-to-order (12-14 weeks)
Target Names	FXN
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.