



CASD1 Antibody

Product Code	CSB-PA846683ZA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96PB1
Immunogen	Recombinant Homo sapiens CASD1 protein
Raised In	Rabbit
Species Reactivity	Homo sapiens
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	N-acetylneuraminate 9-O-acetyltransferase (EC 2.3.1.45) (CAS1 domain-containing protein 1) (Sialate O-acetyltransferase) (SOAT), CASD1, C7orf12
Immunogen Species	Homo sapiens (Human)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	CASD1
Product Modify	Unmodified