



aexT Antibody

Product Code	CSB-PA856653ZA01AXZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q93Q17
Immunogen	Recombinant Aeromonas salmonicida aexT protein
Raised In	Rabbit
Species Reactivity	Aeromonas salmonicida
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	ADP-ribosyltransferase toxin AexT (EC 2.4.2.-) (Exoenzyme T), aexT
Immunogen Species	Aeromonas salmonicida
Lead Time	Made-to-order (12-14 weeks)
Target Names	aexT
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.