



NFYA3 Antibody

Product Code	CSB-PA856683ZA01DOA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q93ZH2
Immunogen	Recombinant Arabidopsis thaliana NFYA3 protein
Raised In	Rabbit
Species Reactivity	Arabidopsis thaliana
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nuclear transcription factor Y subunit A-3 (AtNF-YA-3) (Transcriptional activator HAP2C), NFYA3, HAP2C
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	NFYA3
Product Modify	Unmodified