



PNKP Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA857031ZA01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q96T60 |
| Immunogen | Recombinant Homo sapiens PNKP protein |
| Raised In | Rabbit |
| Species Reactivity | Homo sapiens |
| Tested Applications | ELISA, WB (ensure identification of antigen) |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | Antigen Affinity Purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Bifunctional polynucleotide phosphatase/kinase (DNA 5'-kinase/3'-phosphatase) (Polynucleotide kinase-3'-phosphatase) [Includes: Polynucleotide 3'-phosphatase (EC 3.1.3.32) (2'(3')-polynucleotidase); Polynucleotide 5'-hydroxyl-kinase (EC 2.7.1.78)], PNKP |
| Immunogen Species | Homo sapiens (Human) |
| Lead Time | Made-to-order (12-14 weeks) |
| Gene Names | PNKP |
| Product Modify | Unmodified |