



# CCL19 Antibody

<b>Product Code</b>	CSB-PA859522ZA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q99731
<b>Immunogen</b>	Recombinant Homo sapiens CCL19 protein
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Homo sapiens
<b>Tested Applications</b>	ELISA, WB (ensure identification of antigen)
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	C-C motif chemokine 19 (Beta-chemokine exodus-3) (CK beta-11) (Epstein-Barr virus-induced molecule 1 ligand chemokine) (EBI1 ligand chemokine) (ELC) (Macrophage inflammatory protein 3 beta) (MIP-3-beta) (Small-inducible cytokine A19), CCL19, ELC MIP3B SCYA19
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Lead Time</b>	Made-to-order (12-14 weeks)
<b>Target Names</b>	CCL19
<b>Product Modify</b>	Unmodified
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.