



ERP44 Antibody

Product Code	CSB-PA861125ZA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BS26
Immunogen	Recombinant Homo sapiens ERP44 protein
Raised In	Rabbit
Species Reactivity	Homo sapiens
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Endoplasmic reticulum resident protein 44 (ER protein 44) (ERp44) (Thioredoxin domain-containing protein 4), ERP44, KIAA0573 TXNDC4
Immunogen Species	Homo sapiens (Human)
Lead Time	Made-to-order (12-14 weeks)
Target Names	ERP44
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.