



EXO5 Antibody

Product Code	CSB-PA863995ZA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H790
Immunogen	Recombinant Homo sapiens EXO5 protein
Raised In	Rabbit
Species Reactivity	Homo sapiens
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Exonuclease V (Exo V) (hExo5) (EC 3.1.-.-) (Defects in morphology protein 1 homolog), EXO5, C1orf176 DEM1
Immunogen Species	Homo sapiens (Human)
Lead Time	Made-to-order (12-14 weeks)
Target Names	EXO5
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.