



prpL Antibody

Product Code	CSB-PA875877ZA01EZ X
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HWK6
Immunogen	Recombinant Pseudomonas aeruginosa prpL protein
Raised In	Rabbit
Species Reactivity	Pseudomonas aeruginosa
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Lysyl endopeptidase (EC 3.4.21.50) (Protease IV) (PvdS-regulated endoprotease), prpL
Immunogen Species	Pseudomonas aeruginosa (strain ATCC 15692 / PAO1 / 1C / PRS 101 / LMG 12228)
Lead Time	Made-to-order (12-14 weeks)
Target Names	prpL
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.