



# Lingo1 Antibody

<b>Product Code</b>	CSB-PA880476ZA01MO
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q9D1T0
<b>Immunogen</b>	Recombinant Mus musculus Lingo1 protein
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Mus musculus
<b>Tested Applications</b>	ELISA, WB (ensure identification of antigen)
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Leucine-rich repeat and immunoglobulin-like domain-containing nogo receptor-interacting protein 1 (Leucine-rich repeat neuronal protein 1) (Leucine-rich repeat neuronal protein 6A), Lingo1, Lern1 Lrrn6a
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Lead Time</b>	Made-to-order (12-14 weeks)
<b>Gene Names</b>	Lingo1
<b>Product Modify</b>	Unmodified