



Pburs Antibody

Product Code	CSB-PA893290ZA01DLU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9VJS7
Immunogen	Recombinant Drosophila melanogaster Pburs protein
Raised In	Rabbit
Species Reactivity	Drosophila melanogaster
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Partner of bursicon (Bursicon subunit beta), Pburs, burs-beta
Immunogen Species	Drosophila melanogaster (Fruit fly)
Lead Time	Made-to-order (12-14 weeks)
Target Names	pburs
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.