



# iagB Antibody

|                            |                                                                                  |
|----------------------------|----------------------------------------------------------------------------------|
| <b>Product Code</b>        | CSB-PA315647ZA01SXB                                                              |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                    |
| <b>Uniprot No.</b>         | P0CL15                                                                           |
| <b>Immunogen</b>           | Recombinant Salmonella typhimurium iagB protein                                  |
| <b>Raised In</b>           | Rabbit                                                                           |
| <b>Species Reactivity</b>  | Salmonella typhimurium                                                           |
| <b>Tested Applications</b> | ELISA, WB (ensure identification of antigen)                                     |
| <b>Form</b>                | Liquid                                                                           |
| <b>Conjugate</b>           | Non-conjugated                                                                   |
| <b>Storage Buffer</b>      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| <b>Purification Method</b> | Antigen Affinity Purified                                                        |
| <b>Isotype</b>             | IgG                                                                              |
| <b>Clonality</b>           | Polyclonal                                                                       |
| <b>Alias</b>               | Invasion protein iagB, iagB                                                      |
| <b>Immunogen Species</b>   | Salmonella typhimurium (strain LT2 / SGSC1412 / ATCC 700720)                     |
| <b>Lead Time</b>           | Made-to-order (12-14 weeks)                                                      |
| <b>Target Names</b>        | iagB                                                                             |
| <b>Product Modify</b>      | Unmodified                                                                       |
| <b>Usage</b>               | For Research Use Only. Not for use in diagnostic or therapeutic procedures.      |