




# Omega-hexatoxin-Hv1a Antibody

<b>Product Code</b>	CSB-PA347842ZA01HAD
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P56207
<b>Immunogen</b>	Recombinant Hadronyche versuta Omega-hexatoxin-Hv1a protein (1-37aa )
<b>Species Reactivity</b>	Hadronyche versuta
<b>Tested Applications</b>	ELISA, WB
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	Protein G
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Product Type</b>	Polyclonal Antibody
<b>Immunogen Species</b>	Hadronyche versuta (Blue mountains funnel-web spider) (Atrax versutus)
<b>Research Area</b>	Others
<b>Target Names</b>	N/A

## Image



### Western Blot

Positive WB detected in Recombinant protein  
All lanes: Omega-hexatoxin-Hv1a antibody at 1:2000

### Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution  
Predicted band size: 25 kDa  
Observed band size: 25 kDa