




Omega-hexatoxin-Hv1a Antibody

Product Code	CSB-PA347842ZA01HAD
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P56207
Immunogen	Recombinant Hadronyche versuta Omega-hexatoxin-Hv1a protein (1-37aa)
Species Reactivity	Hadronyche versuta
Tested Applications	ELISA, WB
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein G
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Hadronyche versuta (Blue mountains funnel-web spider) (Atrax versutus)
Research Area	Others
Target Names	N/A

Image



Western Blot

Positive WB detected in Recombinant protein
All lanes: Omega-hexatoxin-Hv1a antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution
Predicted band size: 25 kDa
Observed band size: 25 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.