



Glutenin, low molecular weight subunit 1D1 Antibody

Product Code	CSB-PA319892ZA01TQN
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10386
Immunogen	Recombinant Triticum aestivum Glutenin, low molecular weight subunit 1D1 protein
Raised In	Rabbit
Species Reactivity	Triticum aestivum
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glutenin, low molecular weight subunit 1D1
Immunogen Species	Triticum aestivum (Wheat)
Lead Time	Made-to-order (12-14 weeks)
Product Modify	Unmodified