



# KAE1 Antibody

<b>Product Code</b>	CSB-PA330543ZA01SVG
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P36132
<b>Immunogen</b>	Recombinant Saccharomyces cerevisiae KAE1 protein
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Saccharomyces cerevisiae
<b>Tested Applications</b>	ELISA, WB (ensure identification of antigen)
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	tRNA N6-adenosine threonylcarbamoyltransferase (EC 2.3.1.234) (Kinase-associated endopeptidase 1) (N6-L-threonylcarbamoyladenine synthase) (t(6)A synthase) (t(6)A37 threonylcarbamoyladenine biosynthesis protein KAE1) (tRNA threonylcarbamoyladenine biosynthesis protein KAE1), KAE1
<b>Immunogen Species</b>	Saccharomyces cerevisiae (strain ATCC 204508 / S288c) (Baker's yeast)
<b>Lead Time</b>	Made-to-order (12-14 weeks)
<b>Target Names</b>	KAE1
<b>Product Modify</b>	Unmodified
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.